Øvelse nr. 9: HENFALDSLOVEN

Formål

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- At eftervise henfaldsloven
- At bestemme halveringstiden for $^{137}_{56}\text{Ba}^*$.

Apparatur

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
Lidt ekstra vedr. kilden


General Information

The Cs/Ba-137m Isotope Generator is designed to demonstrate the properties of radioactive decay. Based on the original Union Carbide design, it offers exceptional performance combined with ease of use and safe operation. This product is exempt from USNRC and State licensing and requires no special handling, storage or disposal requirements. A small quantity ($<10\mu\text{Ci}$) of radioactive Cs-137 is


bound on a special ion exchange medium. The Cs-137 parent isotope beta decays with a 30.17y half-life to produce Ba-137m which in turn decays with a 2.55min. half-life, generating a 661.6 keV gamma ray emission. This gamma ray may be readily detected using a GM or scintillation radiation detector. An eluting solution is used to selectively extract the Ba-137 isotope from the exchange medium leaving the parent Cs-137 isotope in place to regenerate more Ba-137. Equilibrium is re-established in less than 1 hour.

Approximately 30 minutes after elution, the residual activity of the Ba-137 solution has decayed to less than 1/1000 of its original activity making it safe for normal disposal.

When used with the eluting solution supplied, bleed through of the Cs-137 parent isotope is less 50 Bq/ml. In order to maintain correct chemical stability, it is important to use only the correct eluting solution. Additional solution may be ordered as part ELSN or prepared by the user as 0.9% NaCl in 0.04M HCl. When making solution, use distilled or DI water to avoid unwanted mineral contamination.

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